

UV HAND HYGIENE TRAINING DEVICE

Helping to eliminate cross contamination

H A N D
I N S P E C T I O N
C A B I N E T S

The Hand Inspection Cabinet will highlight instantly any defects in hand washing techniques and raise awareness of the potential for hygiene problems.



OBJECTIVES :

- To raise awareness in your staff of the importance of good hand hygiene.
- To demonstrate that an effective hand wash requires the correct technique and friction.
- To demonstrate the importance of hands being washed often and well.

HOW IT WORKS :

A fluorescent lotion is applied to the hands and then washed off. When the hands are placed in the cabinet, any remaining lotion will fluoresce under the ultraviolet lamp dramatically illustrating any flaws in the hand washing technique.

DaRo
UV SYSTEMS LTD

T: +44 (0) 1787 881191
E: info@darouv.co.uk



A C C E S S O R I E S



GLITTER BUG PRODUCTS

The Glitter Bug products are formulated from UV fluorescent powder mixed in a special hand lotion so that it does not separate. The pump dispenses just the right amount to avoid wastage, with approximately 200 applications per bottle. Small 20 ml disposable sachets of potion (free trial) are available for demonstration purposes. The potion fluoresces brightly under UV illumination.

Please contact our office for further details and prices of the hand inspection cabinet and its accessories

GLITTER BUG POTION

is used with the Daro Hand Inspection Cabinet

DARO CARRYING CASE

The black and white Carrying Case is specially made to suit the size of the Hand Inspection Cabinet.

Measures: 310 mm Height, 280 mm Depth, 525 mm Width.



PRACTICAL CARRYING CASE

for Daro Hand Inspection Cabinet

ORDER FORM

Hand Inspection Cabinet

Order Code	Description	Price each	Qty.	Sub-total
STMMI	Inspection Cabinet	£163.00		£
STMMI000	Glitterbug Lotion	£16.00		£
STMMI001	Carrying Case	£19.75		£
BLB15	Replacement Lamp	£19.75		£
		Carriage		£
		Sub Total		£
		VAT @17.5%		£
		Total		£

UK Mainland Carriage Costs :- 1-2 cabinets £12.00, 3-5 cabinets £20.00

VAT is chargeable on both goods and carriage

Please call for carriage costs for deliveries outside the UK.

Order values under £50.00 (exclusive of carriage and VAT) will incur an administration charge of £12.00

Organisation Name			
Contact Name			
Telephone Number		Fax Number	
Email Address			
Accounting Address			
Post Code			
Official Order Number			

Please complete this form and send it to DaRo UV Systems, as detailed below.

Churchfield Road, Chilton, Sudbury, Suffolk CO10 2YA

T: 01787 881191 F:01787 374291


E: sales@daroproducts.co.uk W: www.darouv.co.uk

Vat Reg: 898 6614 51 Company Reg. No. 6070435

Bank account number: 50300594 Bank sort code: 20-83-50

HOSPITAL HANDWASHING TECHNIQUE

1



Place solution in palm of wet hand.

2



Massage palm to palm.

3



Rub right palm over left dorsum and vice versa.

4



Rub palm to palm with fingers interlaced.

5



Massage backs of fingers in opposing palm with fingers interlocked.

6



Rotate right thumb clasped in left palm and vice versa.

7



Rotate clasped fingers of left hand in right palm and vice versa.

8



Rinse thoroughly for one minute under running tap water.

DaRo UV Systems

CHURCHFIELD ROAD • CHILTON IND. EST. • SUDBURY • SUFFOLK • CO10 2YA
T: 01787 881191 • F: 01787 374291 • E: sales@daroproducts.co.uk • W: www.darouv.co.uk